


LONGER  
ROTATING SHAFT

MORE REACH  
EASIER CONTROL




**LSOLUTIONS®**

7796 Victor-Mendon Road  
Victor, New York 14564 U.S.A.  
Customer Service: 866.575.3493  
Technical Support: 866.428.9092  
Fax: 585.742.8086  
www.lsisolutions.com

**NEW**  
**COR-KNOT**  
**MINI® DEVICE**

# COR-KNOT MINI® DEVICE

SPINI MINI SHAFT | 17cm Long, Rotating Shaft



Actual Size


LS SOLUTIONS COR-KNOT MINI

SLIDE


SQUEEZE & HOLD

TUG & RELEASE


REMOVE DEVICE



Suture Slot Faces Center



Compression & Tension



Squeeze Lever Only Once








Inspect Fastener



## COR-KNOT MINI® PRODUCT ORDERING

SUPPLIED: STERILE

|  | REORDER    | PRODUCT                       | DESCRIPTION  |
|--|------------|-------------------------------|--|
|  x 6  | REF 031300 | COR-KNOT MINI® DEVICE KIT     | Box of 6 Kits<br>(2 Devices per Kit)                   |
|  x 12   | REF 030850 | COR-KNOT® QUICK LOAD® SINGLES | Box of 12 Singles<br>(1 Fastener per Pouch)            |
|  x 12   | REF 030874 | COR-KNOT® QUICK LOAD® 6-POUCH | Box of 12 Pouches<br>(6 Fasteners per Pouch)           |
|  x 6<br> x 6 | REF 031350 | COR-KNOT MINI® COMBO KIT      | Box of 6 Kits<br>(2 Devices &<br>12 Fasteners per Kit) |

PN031443A 04.20.16

LS SOLUTIONS®