

SEW-EASY[®] DEVICE

Automated Dual Straight Needle Prosthetic Suturing Device

Reliable, effective prosthetic
suture placement with
no exposed needles.
Available in 3.5mm and
5.0mm suture spacing



LS SOLUTIONS[®]

SEW-EASY® DEVICE IN USE












PROSTHETIC SUTURE PLACEMENT

Needles advanced through sewing cuff.



PROSTHETIC SUTURE PLACEMENT

Sutures placed through sewing cuff.

PRODUCT ORDERING			SUPPLIED: STERILE	
	REORDER 3.5mm SPACING	REORDER 5mm SPACING	PRODUCT	DESCRIPTION
 x 6	REF 022200	REF 022797	SEW-EASY® COMBO KIT	Box of 6 Kits (2 Devices, 12 Cassettes, and 2 Valve Handles per Kit)
 x 6	REF 022201	REF 022804	SEW-EASY® DEVICE KIT	Box of 6 Kits (2 Devices and 2 Valve Handles per Kit)
 x 12	REF 022160	REF 022811	SEW-EASY® CASSETTE SINGLES	Box of 12 Cassettes (1 Cassette per Pouch)
 x 6	REF 021900	REF 022580	RAM® DEVICE	Box of 6 Kits (2 Devices per Kit)
 x 12	REF 022039 PLEDGETED REF 022325 NON-PLEDGETED	REF 022668 PLEDGETED REF 022680 NON-PLEDGETED	RAM® COR-SUTURE® QUICK LOAD® SINGLES Pledgeted 2-0 Polyester, Green, Non-Absorbable, 38"	Box of 12 Sutures (1 Suture per Pouch)
 x 12	REF 022047 PLEDGETED REF 022320 NON-PLEDGETED	REF 022671 PLEDGETED REF 022683 NON-PLEDGETED	RAM® COR-SUTURE® QUICK LOAD® SINGLES 2-0 Polyester, White, Non-Absorbable, 38"	Box of 12 Sutures (1 Suture per Pouch)
 x 12	REF 022335 PLEDGETED REF 022340 NON-PLEDGETED	REF 022674 PLEDGETED REF 022686 NON-PLEDGETED	RAM® COR-SUTURE® QUICK LOAD® SINGLES 2-0 Polyester, Green/White Striped, Non-Absorbable, 38"	Box of 12 Sutures (1 Suture per Pouch)
 x 6	REF 022090	REF 022677	RAM® COR-SUTURE® QUICK LOAD® 6 TRAY 2-0 Polyester, White (3) and Green (3), Non-Absorbable, 38"	Box of 6 Trays (6 Sutures per Tray)
 x 12	REF 022540*		SEW-EASY® SNARE	Box of 12 Snares (1 Snare per Pouch)

*Only available in the U.S.

EC REP

EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands



LSI SOLUTIONS®

7796 VICTOR-MENDON RD | VICTOR, NY 14564 USA | CUSTOMER SERVICE: 866.575.3493 | TECH SUPPORT: 866.428.9092 | FAX: 585.742.8086 | WWW.LSISOLUTIONS.COM
The LSI logo, LSI SOLUTIONS, RAM, Sew-Easy, Cor-Suture, and Quick Load are trademarks and registered trademarks of LSI Solutions, Inc. Copyright © 2017, LSI SOLUTIONS®. All Rights Reserved.

PN021565D 04-21-23